

## Standard Operating Procedure for Hazardous Chemicals

**Principal Investigators:** Chung-Jui Tsai and Scott A. Harding  
**Building and rooms:** Davison Life Sciences Building, Lab B310

<b>Chemical(s)</b>	Tris-HCl, EDTA (ethylenediaminetetraacetic acid), $\beta$ ME ( $\beta$ -mercaptoethanol), SDS (Sodium Dodecyl Sulfate), Ethanol, KOAc (potassium acetate), Isopropanol, liquid nitrogen (LN), sodium chloride (NaCl).
<b>Process</b>	<b>Plant Genomic DNA extraction - miniprep</b>
<b>Specific Hazards</b> <i>referred to MSDSs for more detailed information</i>	<u>Isopropanol</u> , <u>Ethanol</u> : Flammable. <u><math>\beta</math>ME</u> : flammable, and may be harmful upon skin/eye contact, ingestion or inhalation. <u>LN</u> : frostbite hazard. <u>KOAc</u> : Hazardous in case of eye contact (irritant).
<b>Personal protective equipment</b>	<input checked="" type="checkbox"/> 3-5 mil nitrile gloves <input checked="" type="checkbox"/> double gloves (w/ concentrated stock) <input checked="" type="checkbox"/> lab coat (except when used in microcentrifuge tubes) <input checked="" type="checkbox"/> chemical safety goggles (when splash potential exists)
<b>Engineering/ventilation controls</b>	All operations involving <u><math>\beta</math>ME</u> must be done in a chemical fume hood.
<b>Special handling procedures and storage requirements</b>	Store <u>isopropanol</u> , <u><math>\beta</math>ME</u> in the flammable cabinets under the hood B. LN: Store and use with adequate ventilation. Under normal conditions these containers will periodically vent product. Do not plug, remove, or tamper with pressure relief device.
<b>Spill and accident procedures</b>	<u>Skin exposure</u> : Rinse affected skin with plenty of water while removing contaminated clothing/shoes. Rinse for > 15 minutes.  <u>Eye exposure</u> : Wash eyes for > 15 minutes. For both cases, seek medical attention immediately.
<i>for hazardous chemicals only</i>	<u>Small</u> (< 2L): Absorb with vermiculite or spill pads and transfer absorbed material to a closed container. Label and date as hazardous waste for disposal. Notify PIs.  <u>Large</u> (> 2L): Evacuate the room, notify PIs and call 2-5801 to request emergency spill assistance from the Environmental Safety Division.
<b>Waste disposal</b>	Dispose waste as regular disposing method.
<b>Special approval</b>	No special authorization needed after SOP training & reading MSDSs.
<b>Prepared by</b>	Name/date: Kate Tay, 9/2/2009; updated by CJ Tsai 6/27/2013
<b>Reviewed by</b>	Name/date: C-J Tsai, 9/4/2009; 6/27/2013

## Plant Genomic DNA (gDNA) extraction miniprep v2.2

### Reagents and Buffers:

<u>gDNA Extraction Buffer:</u>	[Stock]	[Final]	For 50 ml
Tris-HCl (pH 8)	1.0 M	50 mM	2.5 ml
EDTA (pH 8)	0.5 M	10 mM	1.0 ml
NaCl	5.0 M	100 mM	1.0 ml
SDS	20%	1%	2.5 ml
ddH <sub>2</sub> O	--	--	43.0 ml
β-mercaptoethanol	--	10 mM	34 μl (add right before use)

### Potassium acetate (KOAc) pH4.8:

	For 50 ml
ddH <sub>2</sub> O	15 ml
Acetic acid glacial	15 ml * <b>Always add acid into water, NOT water into acid</b>
KOH pellet	to bring pH up to pH = 4.8, then bring volume to 50 ml.

\* **Maintain the bottle in a cold water bath to dissipate the heat during the prep and pH measurement.**

**Equipment:** 65°C water bath

### Cleaning Micro-pestle:

Sonicate in 70% EtOH for 20 min, rinse with plenty of ddH<sub>2</sub>O, dry and (optional) autoclaved.

### **Procedures:**

1. Collect tissue (e.g., a single punch, a small apical leaf or ~ 50 mg) into a 1.5 ml microcentrifuge tube and keep on ice if multiple samples are to be collected. Snap-freeze in LN immediately after return to the lab. Store at -80°C until use.
2. Grind tissue to a fine powder using a micro-pestle under liquid nitrogen (fill to one-third of the tube). Proceed to the extraction step immediately or store samples at -80 °C until use.
3. Add 750 μl gDNA extraction buffer to each tube of frozen powder. Mix well by vortexing.
4. Incubate @ 65°C for 20 min. Invert the tubes to mix every 5 min.
5. Add 200 μl **ice-cold** 5M KOAc and Invert 5 times to mix. Incubate on ice for 20 min.
6. Centrifuge at max speed for 10 min at room temperature.
7. While waiting, add 480 μl (or 0.6 Vol. of DNA extract) isopropanol into new 1.5ml tubes.
8. Transfer 800 μl of the extract to the tube containing isopropanol. Invert 5 times to mix.
9. Centrifuge at max speed for 10 min at room temperature to pellet DNA.
10. Wash with 500 μl of 70% EtOH, invert the tube several times to mix, and centrifuge for 5 min. Repeat the wash once more.
11. Vacuum dry using centrivap for ~10min.
12. Resuspend DNA in 50 μl ddH<sub>2</sub>O. If the level of co-purifying RNA is high, resuspend DNA in ddH<sub>2</sub>O containing 10 μg/ml RNase A (make a master mix, just enough for your samples) and incubate at 37°C for 20 min.
13. QC 1 μl on a 1% agarose gel. Some amount of RNA usually co-purify with the gDNA, therefore, QC by Nanodrop is not necessary/meaningful (unless RNaseA treatment & EtOH ppt are done). Estimate gDNA concentration based on the ladders. You need 1-5 ng for each PCR reaction. If yield is high, dilute gDNA 10X and check on a gel again.

### **References:**

Stephen L. Dellaporta, Jonathan Wood, James B. Hicks. (1983) A plant DNA miniprep: Version II. Plant Molecular Biology Reporter, Vol.1, 4:19-21