

Standard Operating Procedure for Hazardous Chemicals

Principal Investigators: Chung-Jui Tsai and Scott A. Harding
Building and rooms: Davison Life Sciences Building, Lab B310

Chemical(s)	Potassium ferrocyanide, Potassium ferricyanide, EDTA (ethylenediaminetetraacetic acid), Ascorbic Acid, PVP (Polyvinylpyrrolidone), X-Gluc (5-bromo-4-chloro-3-indolyl-beta-D-glucuronic acid), DMF (Dimethylformamide), Acetone, Ethanol.
Process	GUS Staining Protocol
Specific Hazards referred to MSDSs for more detailed information	<u>Ascorbic Acid</u> : Slightly hazardous. Irritant. <u>PVP</u> : Irritant. Hygroscopic. Do not breathe dust. <u>DMF</u> : Irritant. Toxic. Possible reproductive toxin. <u>Acetone, Ethanol</u> : Flammable. X-Gluc: Irritant.
Personal protective equipment	<input checked="" type="checkbox"/> 3-5 mil nitrile gloves <input checked="" type="checkbox"/> double gloves (w/ concentrated stock) <input checked="" type="checkbox"/> lab coat (except when used in microcentrifuge tubes) <input checked="" type="checkbox"/> chemical safety goggles (when splash potential exists)
Engineering/ventilation controls	Emergency shower and eyewash accessible. Chemical fume hood (when working with large quantities)
Special handling procedures and storage requirements	<u>X-Gluc</u> : stored @ -20°C, light sensitive. <u>K-ferrocyanide, K-ferricyanide</u> : stored @ 4°C, light sensitive. <u>PVP</u> : stored @ 4°C. <u>DMF, Acetone, Ethanol</u> : flammable cabinet, Rm B310, store away from oxidizers. <u>Ascorbic Acid</u> : Air and light sensitive.
Spill and accident procedures	<u>Skin exposure</u> : Rinse affected skin with plenty of water while removing contaminated clothing and shoes. Rinse for at least 15 minutes. Seek medical attention, as appropriate. <u>Eye exposure</u> : Wash eyes for at least 15 minutes, lifting the upper and lower eyelids. Seek medical attention, as appropriate.
<i>for hazardous chemicals only</i>	<u>Small</u> (< 2L): Absorb with vermiculite or spill pads and transfer to a closed container. Label and date as hazardous waste for disposal. Notify PIs. <u>Large</u> (> 2L): Evacuate the room, notify PIs and call 2-5801 to request emergency assistance from Environmental Safety Division
Waste disposal	<u>DMF</u> , and <u>Acetone</u> wastes must be collected and labeled as hazardous waste according to the SOP for Hazardous Waste Disposal.
Special approval	No special authorization needed after SOP training & reading MSDSs.
Prepared by	Name/date: C-J Tsai, 8/8/09
Reviewed by	Name/date: Kate Tay, 9/2/09

GUS Staining Protocol

Reagents and Buffers:

1. 50 mM Potassium ferrocyanide (M.W. = 422.39): 0.2110 g in 10 mL of ddH₂O, Stored @ 4°C, Dark.
2. 50 mM Potassium ferricyanide (M.W. = 329.25): 0.1645 g in 10 mL of ddH₂O, Stored @ 4°C, Dark.
3. 1 M Potassium phosphate buffer (pH 7):
3.9 mL of 1 M monobasic (KH₂PO₄) + 6.1 mL of dibasic (K₂HPO₄). **Don't use a pH meter.**
4. 0.5 M EDTA
5. 0.5M Ascorbic Acid – Sodium salt (M.W. = 198.1):
0.9905g in 10ml of ddH₂O, aliquot in amber tube, Stored @ -20°C, Dark.
6. 10% PVP: 1g in 10 ml of ddH₂O, Stored @ 4°C, Dark.
* Use ONLY the soluble PVP (Sigma PVP40T), NOT the poly-PVP.
7. 20 mM X-Gluc (M.W. = 521) in DMF:
Add 1.92 mL of DMF directly to the 100 mg vial of X-gluc vial to make 100 mM stock.
Dilute 5X with DMF to make 20 mM stock. Stored @ -20°C, Dark.
X-Gluc is expensive; do not attempt to weight a sub-portion from the vial (wasteful).
8. Triton X-100
9. 90% Acetone (dilute before use)
10. 70%, 50%, 30%, EtOH.

X-Gluc Staining buffer: **FRESH** (1ml for 4x 96-well)

<u>Component</u>	<u>Stock</u>	<u>Final</u>	<u>For 1 mL</u>
ddH ₂ O	--	--	759 µL
K-phosphate, pH7	1 M	100 mM	100 µL
EDTA	0.5 M	10 mM	20 µL
Triton X-100	--	0.1%	1 µL (pipet slowly)
PVP #	10%	2%	200 µL
K-ferrocyanide [#]	50 mM	0.5 mM	10 µL
K-ferricyanide [#]	50 mM	0.5 mM	10 µL
Ascorbic Acid – Na salt ^{##}	0.5M	25mM	50 µL
X-Gluc ^{##}	20mM	1mM	50µL

Stored @ 4°C, Dark,

Stored @ -20°C, Dark.

Washing buffer: 100 mM K-phosphate, pH 7 (dilute 1 mL of 1M stock with 9 mL of ddH₂O)

Procedures:

1. Use a sharp razor to dissect tissues into pieces of less than 5 mm in any dimension.
2. Submerge immediately in 90% acetone (to fix protein) for **exactly** 30 min (use a timer to record the interval between each tissue dissection when processing multiple samples).
3. Rinse thoroughly with 100 mM potassium phosphate buffer (pH7), 3 times.
4. Add enough (e.g., 250µl) of **FRESH** X-Gluc staining buffer into an appropriate microplate (e.g., 96-well format) and make sure tissues are submerged in the staining solution (while conserving X-gluc).
5. Apply a mild/water vacuum for 10 min to facilitate infiltration.
6. Incubate @ 37°C for 12 hrs (monitor periodically, stop with sign of tissue browning)

7. Discard staining solution and incubate tissues in 70% EtOH to remove chlorophyll at room temp or 37°C. Change the EtOH solution once or twice until all chlorophyll is clear. Tissues can be stored in 70% EtOH (easy to dry out) or water (may get moldy).
8. Before doing vibratome sectioning, rehydrate the tissues through 50% EtOH → 30% EtOH → ddH₂O gradually, 1 hr each (to avoid damaging tissue structure).
9. Make cross section of 70-100 μm thick using a vibratome and seal the slide with nail polish
10. Observe under a microscope.