

Standard Operating Procedure for Hazardous Chemicals

Principal Investigators: Chung-Jui Tsai and Scott A. Harding

Building and rooms: Davison Life Sciences B310

Chemical(s)	1-Butanol, Hydrochloric Acid, Chloroform, Ferric Ammonium Sulfate, Methanol
Process	Condensed Tannin Assay
Specific Hazards <i>referred to MSDSs for more detailed information</i>	<u>Hydrochloric acid</u> is corrosive and liquid or vapor may cause burns to skin/eye. <u>1-Butanol</u> is flammable, and may be harmful upon skin/eye contact, ingestion or inhalation. <u>Chloroform</u> is a toxic solvent and contact to skin or eyes is dangerous. <u>Butanol</u> fumes are toxic. <u>Methanol</u> is flammable and harmful upon ingestion
Personal protective equipment	Must wear 3-5 mil nitrile gloves. Chemical safety goggles and lab coat should be worn when splash potential exist
Engineering/ventilation controls	All operations involving chloroform, hydrochloric acid and butanol must be done in a chemical fume hood.
Special handling procedures and storage requirements	Store <u>chloroform</u> , <u>1-butanol</u> and <u>methanol</u> in the flammable cabinet, and <u>hydrochloric acid</u> in the designated cabinet underneath the chemical hood in Rm <u>B310</u> . Store all the other buffers on bench top.
Spill and accident procedures <i>for hazardous chemicals only</i>	<p><u>Skin exposure</u>: Rinse affected skin with plenty of water while removing contaminated clothing/shoes. Rinse for > 15 minutes.</p> <p><u>Eye exposure</u>: Wash eyes for > 15 minutes.</p> <p>For both cases, seek medical attention immediately.</p> <p><u>Small</u> (< 2L): Absorb with vermiculite or spill pads and transfer absorbed material to a closed container. Label and date as hazardous waste for disposal. Notify PIs.</p> <p><u>Large</u> (> 2L): Evacuate the room, notify PIs and call 2-5801 to request emergency spill assistance from Environmental Safety Division.</p>
Waste disposal	Chloroform and 1-butanol waste must be collected and labeled as separate hazardous wastes according to the SOP for Hazardous Waste Disposal.
Special approval	No special authorization needed after SOP training & reading MSDSs.
Prepared by	Name/date: SA Harding, 8/18/2008.
Reviewed by	Name/date:

Condensed Tannin (CT) Assay Method

Solutions and Reagents needed

Butanol-HCl mix (5% HCl in 1-butanol): Slowly add 10 ml concentrated hydrochloric acid (HCl) to 190 ml 1-butanol in a glass-stoppered storage bottle in the fume hood. Swirl to mix.

FAS (2% ferric ammonium sulfate in 2M HCl): In the fume hood, add 2 ml HCl to 10.1 ml H₂O. Mix, and then add 210 mg ferric ammonium sulfate until completely dissolved.

BFAS (Butanol-FAS mix): Just prior to use, in the fume hood, add 1 ml of FAS to 50 ml of the Butanol-HCl mix. Mix well.

Chloroform, Methanol, heating block

Precautions: Gloves should be worn throughout the procedure. Goggles should be worn when working with chloroform, butanol or hydrochloric acid.

Procedure:

1. Leaf discs (fresh or freeze-dried) or tissue powder (freeze-dried) can be used. Use a standard paper hole puncher (1/4") to collect leaf discs into an Eppendorf tube containing methanol (100ml/disc). 2-4 discs, or 2-10 mg freeze-dried powder should be sufficient. Sonicate in a bath sonicator for 30 min. Transfer the de-pigmented discs to a new tube. Re-extract the tissue if green remains.

If the tissue being measured does not have chlorophyll (e.g., roots or twigs), the following chlorophyll removal steps (#2) can be skipped.

2. Add 1.5 Vol chloroform and 2 Vol H₂O (**chloroform first, then H₂O**) to the methanol extract. Vortex and centrifuge at 10-15,000xg for 1 min to separate phases. Transfer the clear supernatant to a new tube, and dispose the lower chloroform phase as organic waste. Add small amounts of H₂O if separate phases do not form.
3. Evaporate the de-pigmented discs/powder (containing "insoluble" CTs) and the clear supernatant fraction (containing "soluble" CTs) to dryness using a Centrivap or equivalent.
4. CT color assay. To dried tissue and/or extract (two different assays), add 500ml of the BFAS solution. Close the microtube lid and attach a lid-lock. Heat the tube at 95°C for 20 min. Cool to room temperature and read Abs₅₅₀. If quantitative rather than comparative results are desired, a standard curve using purified CTs will have to be prepared. If the color is dark, dilute accordingly using Butanol-HCl (NOT BFAS, water or methanol) and repeat the measurement.